

# Sizing Chart for Sigvaris Sheer Fashion 15-20 mmHg Level of Compression



## Sigvaris 120 Sheer Fashion Maternity Pantyhose

Use pre-pregnancy weight and match your height and weight to find the sizing range you fall into.

lbs	95	100	105	110	115	120	125	130	135	140	145	150	155	160	165	170	175	180	185	190	195	200	205	210	215	220	225	230	235	240	245	250	
Kg	43	45	48	50	52	54	57	59	61	64	66	68	70	73	75	77	80	82	84	86	89	91	93	95	98	100	102	105	107	109	111	114	
4'10" 147cm																																	
4'11" 150cm																																	
5' 152cm																																	
5'1" 155cm																																	
5'2" 157cm																																	
5'3" 160cm																																	
5'4" 162cm																																	
5'5" 165cm																																	
5'6" 167cm																																	
5'7" 170cm																																	
5'8" 173cm																																	
5'9" 175cm																																	
5'10" 178cm																																	
5'11" 180cm																																	
6' 183cm																																	

15-20 mmHg is applied to an average ankle circumference for sizes:

**A & B:** 7.5" - 8.5" (19 - 22cm)

**C & D:** 8.5" - 10" (22 - 25cm)

**E:** 10" - 11.5" (25 - 29cm)

*only at LegSmart.com*